

Manage Test Panel for Tests and Results


Making tests and results entry for specified lab easier

[Contact Species360 Support](#)


In tests and results in ZIMS for Medical, you can provide a list of known laboratory panels that are to be used for a specific sample type.

To save data entry time you can use Test Panels by selecting Manage Favorite Test Panel. Commonly used Test Panels will be preloaded to select from.


Test panels can be added from the Test and Results entry screen from the Manage Favorite Test Panel:



Additionally, this can be managed from the start menu, Tools, and Manage Test Panels.




Once you have loaded the test panel management screen, you can create your new Test Panels by selecting Add New Panel.




The Panel Name must be unique within your facility. The Valid Sample Type is a multiple select dropdown that allows you to select more than one type.

You are also able to select an external laboratory or your institution to run the test. If the external lab you are selecting is not already in the global ZIMS institution list, you must add it here.

It is the mission of Species360 to facilitate international collaboration in the collection and sharing of information on animals and their environments for zoos, aquariums and related organizations.



In the Manage Test Panel grid, hovering over the right hand side pop out will display what tests you have created for that Panel.



Laboratory	Panel Name	Sample type(s)
Bronson Animal Disease Diagnostic La...	test chem	Serum,Plasma,Whole Blood
Charles River UK Ltd	Serum 1	Serum
Audubon Zoo	zebra blood	Whole E
Bronson Animal Disease Diagnostic La...	pn1	Serum

Test Item: Hematocrit (manual), White Blood Cell count (Natt-Herrick.manual), Lymphocyte absolute count (manual), Basophil absolute count (manual), Partial Pressure of Oxygen

To use the Test Panel select the appropriate sample and complete the mandatory fields in the Test Request Details. Select the desired Test Panel from the dropdown. Then select Insert Below.

Add Test Result

➕ Add New Clinical Diagnosis

Test Request Details

Date Requested * Mar 28, 2014 **Requested By *** Adrienne Miller

External Laboratory My Institution

Analysis By Adrienne Miller **Analysis Start Date & Time** Mar 28, 2014 00:00

Analysis Equipment

Insufficient sample

Sample Info

Sample Type Whole Blood

Sample Collection Date & Time Mar 01, 2014 00:00

Additives/Preservatives

Anatomical Source/Tissue

Sample Quality Edit

Color

Color Intensity

Clarity

Consistency

Additional Characteristics clots present

Degraded No

Test Results

Test Panel -- Select One -- Insert below ⚙

Notes/Comments

1. zebra blood

The tests you have recorded for your Panel will display in the order that you arranged them in the Test Panel you created.

➕ Add New Clinical Diagnosis

Test Request Details

Date Requested * Mar 28, 2014 **Requested By *** Adrienne Miller

External Laboratory My Institution

Analysis By Adrienne Miller **Analysis Start Date & Time** Mar 28, 2014 00:00

Analysis Equipment

Insufficient sample

Sample Info

Sample Type Whole Blood

Sample Collection Date & Time Mar 01, 2014 00:00

Additives/Preservatives

Anatomical Source/Tissue

Sample Quality Edit

Color

Color Intensity

Clarity

Consistency

Additional Characteristics clots present

Degraded No

Test Results

Test Panel zebra blood Insert below ⚙

Add New Test

Test	Primary Result	Secondary Result	Low RI	High RI	Evaluation	Excl'd. from RI	Clinical Finding	Reviewed	Remove
HCT	% (c)	N/A	~	~	~			<input type="checkbox"/>	✖
Natt-Herrick WBC	*10 ³ cells/μL	N/A	~	~	~			<input type="checkbox"/>	✖
Lymphocytes	*10 ⁹ cells/L	N/A	~	~	~			<input type="checkbox"/>	✖
Basophils	*10 ⁹ cells/L	N/A	~	~	~			<input type="checkbox"/>	✖
pO2	mmHg	N/A	~	~	~			<input type="checkbox"/>	✖

Notes/Comments

Now record the test results. If the default Unit of Measure is not what you desire it to be you can use the dropdown selections to select the appropriate unit.

It is the mission of Species360 to facilitate international collaboration in the collection and sharing of information on animals and their environments for zoos, aquariums and related organizations.

Test Results	
Test Panel	zebra blood
Test	Primary Result
HCT	% (c)
Natt-Herrick WBC	*10 ³ cells/μL
Lymphocytes	2.64 *10 ⁹ cells/L
Basophils	76 -- Select One --
pO2	86.5 *10 ⁹ cells per litre
<input type="checkbox"/> Notes/Comments	

Revised 10 August 2017

** Species360 Organizational name change added on 07/18/2016*

It is the mission of Species360 to facilitate international collaboration in the collection and sharing of information on animals and their environments for zoos, aquariums and related organizations.